

## Cs-corrected High-Resolution Analytical STEM


T. Oikawa<sup>1</sup> and M. Kawasaki<sup>2</sup>

<sup>1</sup> JEOL(Europe) SAS, Espace Claude Monet-1, Allée de Giverny, 78290 Croissy-sur-Seine, France ; [oikawa@jeol.fr](mailto:oikawa@jeol.fr)

<sup>2</sup> JEOL USA, Inc., 11 Dearborn Road, Peabody, MA 01960, USA

Recently, scanning transmission electron microscopy (STEM) has made huge progress owing to a stable field emission gun (FEG) [1] and sophisticated lens designing technologies [2]. Also, a spherical aberration corrected (Cs-corrected) lens [3] is now becoming widely spread as commercial instrument. In this paper, very high-resolution (atomic scale) elemental analysis was successively applied to a precipitate in the metal, with Cs-corrected STEM.

An instrument of Atomic Scale Analytical Electron Microscope is the JEM-2100F (200 kV FEG TEM/STEM) equipped Cs-corrector for probe forming lens. The data introduced in figure 1 was obtained from Si<sub>3</sub>N<sub>4</sub> ceramics using this microscope. Figure 1(a) shows a high-resolution high angle annular dark-field (HAADF) image obtained with Cs-corrected STEM in the JEM-2100F. The data was obtained at beam diameter of 0.1 nm, and it shows high quality of S/N ratio owing to a high beam intensity of the Cs-corrected STEM. Figure 1(b) shows elemental maps of Si-L and N-K, reconstructed from data series of EELS spectrum-imaging for corresponding energies. Owing to fine beam and higher beam intensity of the Cs-corrected STEM, the elemental signals were well detected and the maps show atomic order resolution. Figure 1(c) shows atomic location model of Si<sub>3</sub>N<sub>4</sub> ceramics projected at direction for incident electron beam orientation for a reference. It is ascertained that the elemental maps are similar to the atomic model.



**Figure 1.** Atomic scale resolution of elemental mapping with EELS spectrum-imaging in Cs-corrected STEM. A specimen was Si<sub>3</sub>N<sub>4</sub> ceramics. (a) A HAADF image of Si<sub>3</sub>N<sub>4</sub> ceramics; (b) Elemental maps of Si and N constructed by EELS spectrum-imaging, showing with color overlapped image and magnified HAADF image; (c) Atomic location model projected at incident beam orientation.

Atomic order resolution EELS-maps of Si<sub>3</sub>N<sub>4</sub> ceramics were reconstructed from a data series obtained by the EELS spectrum-imaging (SI) method. It is ascertained that the JEM-2100F / Cs-corrector is able to realize an atomic resolution elemental mapping.

[1] T. Honda *et al.*, Ultramicroscopy, 54, (1994), 132.

[2] K. Tsuno, J. Electron microscopy, 48, (1999), 801.

[3] M. Haider *et al.*, Nature, 392, (1998) 768.